Continuous Wave II Arthroscopy Pump

Operating Instruction Manual

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Mode d'emploi

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CONTINUOUS WAVE II ARTHROSCOPY PUMP

Operating Manual
Continuous Wave II Arthroscopy Pump

English	1
Deutsch	17
Français	33
Italiano	49
Español	65

Operating Manual
Continuous Wave II Arthroscopy Pump

AR-6450
CE 0124

GB

This operating manual provides important information for the safe use of all components of the Arthrex Continuous Wave II Arthroscopy Pump AR-6450 including accessories. Read the manual thoroughly prior to using this system and keep it in an easily accessible place.

Follow all safety precautions!

The Continuous Wave II Arthroscopy Pump is designed and tested in accordance with DIN EN 60601-1.

According to DIN EN 60601 this product is type BF, class 1.

According to MDD93/42/EEC, Annex IX, Rule 11, this product is classified as a class IIa product.

Prior to use in a surgical procedure, carefully unpack and inspect the components for any signs of damage, which may have occurred during shipping. If shipping damage is suspected, notify Arthrex or the authorized Arthrex distributor immediately. Any damage may compromise patient safety.

If any transport or first installation damage is not reported in writing within 7 days of receiving the product, the warranty could be rendered void.

Please refer to our general terms of business.

CONTENTS	PAGE
1 INTRODUCTION	4
2 PRODUCT DESCRIPTION	5
3 SET-UP AND OPERATION	7
3.1 Set-up	7
3.2 Connecting the tubing	8
3.3 Operation	9
3.4 Using the remote control	10
3.5 Safety features of the pump	11
4 CLEANING AND DESINFECTION	12
5 MAINTENANCE	13
5.1 Troubleshooting	13
5.2 Replacing the fuses	14
5.3 Guarantee	14
6 TECHNICAL DATA	15

Special instructions are emphasized as follows :

NOTE: *This term contains important information regarding set-up and operation to facilitate use and obtain effective results.*

This term contains critical information regarding safe handling and use of this system. Device malfunction or property damage could result if these instructions are not followed.

WARNING: *This term contains critical information by identifying conditions or practices that may result in injury or loss of life, if these instructions are not followed.*

1 INTRODUCTION

The Arthrex Continuous Wave II Arthroscopy Pump AR-6450 is a safe, reliable, user-friendly system which maintains constant, non-pulsed control of intra-articular rinsing and distension pressure throughout all phases of an arthroscopic surgical procedure.

The integrated pressure sensor can be directly connected to the inflow cannulas for direct pressure measurement. There is no need for complex secondary measuring systems or extra venous access.

The continuous pulse-free flow reacts immediately to changes in the intra-articular pressure so that joint distension can be sustained even under high shaver extraction volumes or secondary outflow.

The nominal values for flow and pressure are adjusted at the control unit keyboard or by using the remote control. The unit has two different tubes (a. roll pump tube, b. patient tubing), so that the roll pump tubing can be economically used for an entire day of surgery

WARNING *The patient tubing must be replaced before each new surgical procedure.*


The control unit is designed for continuous operation. The system is protected by a resetting thermo fuse. In the event of a fault, e.g. excess current, the fuse heats up suddenly and the device switches off. As soon as the thermoelement cools down, the fuse resets automatically in about 30 seconds. The device can be switched on and started again.

WARNING *This device should only be used by trained, qualified medical staff. Use of this equipment by untrained personnel or use of this equipment for indications other than those described may cause patient injury.*

2 PRODUCT DESCRIPTION


Control unit – front view

The front panel of the control unit AR-6450 has the following elements and symbols:

Control unit – rear panel

The rear panel of the control unit has the following functions and symbols:

3 SET-UP AND OPERATION

3.1 Set-up

Place the control unit on a flat, dry surface, such as an arthroscopy equipment cart or a small instrument table.

Connect the main cable lead with the main jack on the back panel of the device. Connect the cable plug to the power socket. If an extension cord is required, it must be grounded. Once the control unit has been connected, it is ready for use. The roll pump is in stand-by mode (the green LED in the stand-by button on the front panel flashes). Simply press the standby button to turn the roll pump on (the green LED in the standby button on the front panel stays on). Press standby again to return it to stand-by mode.

NOTE: *The system will run continuously without a fault from power failures for up to 50 minutes. If a longer power failure occurs, the system will return to its basic mode.*

If high-frequency instruments are in use or defibrillation of the patient is required, ensure equipment is not in direct contact with the patient.

3.2 Connecting the tubing

The roll pump tubing should be connected by the circulator.

- To connect the tubing, release the locking device in the roll pump housing and open it completely. Swivel the cover plate upwards: the roller mechanism is now free.
 - Place the elastic part of the roll pump tubing on the side near the rinsing fluid input into the upper device of the tubing guide and guide the elastic part of the tubing over the rollers. Place the output side of the elastic section into the lower area of the tubing guide, close the cover plate and lock.
- The tubing is connected properly when the green clip ring of the tubing is aligned with the green mark on the guide.
- Close and lock the cover plate.

NOTE: *The locking device must lock securely. A switch in the cover plate prevents the machine from operating if the mechanical section is not locked securely.*

- Connect pressure sensor to the jack.
- Puncture the rinsing fluid bags with the needles on the tubing. If only one rinsing fluid bag is being used, seal the second rinsing fluid line with an additional clip.
- Then turn the ON/standby switch to ON.
- The surgical staff takes the sterile patient tubing out of the sterile pack and hands the connection for the BaktoSafe connector for the roll pump tubing to the circulator.
- The circulator connects the two tubing systems.

WARNING: *The patient tubing must be changed for each patient.*

NOTE: *Every roll pump tubing is supplied with a sterile BaktoSafe cap. The roll pump tubing can be used throughout the entire day of surgery.*

3.3 Operation

The following pressure levels are recommended for surgery:

- Knee arthroscopy 35 – 60 mmHg
- Shoulder arthroscopy 60 – 80 mmHg
- Arthroscopy in small joints 50 – 70 mmHg

The ideal intra-articular pressure depends on the indications for the arthroscopic procedure, bleeding tendency and the possibility of ischemia.

WARNING: *Although the safety and effectiveness of the pumps is verified and documented, the risk of extra-articular edemas must always be kept in mind for patients with pathologically changed articular capsules and for procedures involving an opening of the capsule (e.g. lateral release). Always start with the lowest possible pressure to achieve the desired joint distension. Continue to increase distension pressure until a clear liquid medium is obtained.*

1. After adjusting the required pressure, remove the cap from the distal end of the patient tubing.
2. Open all tubing clips.
3. Switch the pump on.
4. Fill the entire length of the hose to remove any bubbles

*After the air has been purged out of the system, close the clip at the sterile end of the system. The pump must turn off when the tubing is clamped to ensure that the sensor system is working properly.
If the pump does not turn off, check that the clip is firmly closed.
If the pump continues to run and the sensor does not work properly, replace the roll pump tubing.*

5. Connect the patient tubing to the inflow cannula.

NOTE: *The high flow arthroscopy sleeve AR-3031 should be used for optimum flow when rinsing through the inflow cannula.*

6. Open the clip on the patient tubing to release the flow.
7. Once the set pressure is reached, the pump switches off. When the pressure drops, the pump automatically starts again and pumps rinsing fluid into the joint at the set pressure and flow.
8. When the procedure is completed, close all clips and turn the pump off.

3.4 Using the remote control

The Arthrex Continuous Wave II Arthroscopy Pump can also be operated using the remote control (AR-6306) so that the surgeon can adjust flow rate and pressure directly from his operating position.

*The remote control **CANNOT** be sterilized. The remote control and its cable must be covered with a sterile cover (part number AR-6307).*

To set up the remote control, put the plug into the jack on the back of the pump. Cover the device with the sterile cover, ensuring that the full length of the cable is covered.

***ALWAYS** hold the plug to disconnect the cable, otherwise the cable can be damaged.*

3.5 Safety features of the pump

Abnormal operation

The pump has internal circuits to ensure the pressure sensor circuit remains within normal limits. In the event of a fault, the pump will switch off and the decimal point in the three-digit pressure display will flash.

If abnormal operation occurs, disinfect the pump, pack it in the original packaging and send it to Arthrex, accompanied by a brief description of the malfunction.

Overpressure sensor

The overpressure sensor circuit in roll pump AR-6450 measures the pressure of the liquid in the system. Overpressure can occur when the flow is abruptly interrupted or the joint is suddenly positioned in a way which reduces the joint capsule (e.g. bending the knee joint in the "figure 4" position). If overpressure occurs, a red control lamp lights up on the front of the pump and an intermittent acoustic alarm sounds. Pump operation is interrupted until the pressure returns to the set range. Pressure in the joint can be reduced by opening an outflow and by moving the joint to another stress-free position.

Locking the roll pump housing

The pump switches off when the roll pump housing is opened. The locking mechanism prevents any access to the rotating parts while the device is operating.

Fuse

The fuse protects the pump in the event of a fault. In the case of excess current or a short-circuit, the fuse heats up and the pump switches off automatically. The fuse cools down again after approximately 30 seconds. The pump can then be switched on and started up again.

If the pump won't switch on again even after waiting for a longer period of time, disinfect the pump, pack it in the original packaging and send it to Arthrex, accompanied by a brief description of the malfunction.

4 CLEANING AND DESINFEKTION

Control unit

The control unit AR-6450 can be cleaned and disinfected using commercially available surface disinfectants listed.

Always comply with the instructions issued by the manufacturer of the disinfectant.

NEVER use liquid to clean the connections on the remote control and main cable. Remove dust regularly with dry compressed air.

Remote control

The remote control AR-6450 can be cleaned and disinfected using commercially available surface disinfectants listed.

Always comply with the instructions issued by the manufacturer of the disinfectant.

NEVER use liquid on the remote control or other components. Liquid on the connector pins can cause faults in the equipment.

5 MAINTENANCE

Except for replacing the fuses, the Arthrex Continuous Wave II Arthroscopy Pump needs no special maintenance. If the roll pump should malfunction, please contact Arthrex.

5.1 Troubleshooting

Arthrex customer service should be informed immediately in the event of any damage or malfunctions in the equipment.

Fault	Causes and remedies
No power	<ol style="list-style-type: none"> 1. Check the main cable. 2. Send in for repairs.
The pump doesn't work when it is switched on.	<ol style="list-style-type: none"> 1. Check that the housing is closed and locked. 2. Check if the tubing between the cover flap and rollers is jammed. 3. Check that the pressure sensor is connected correctly. 4. Check whether the tubing clips and shut-off valves are open.
No or inadequate flow	<ol style="list-style-type: none"> 1. Check that all tubing clips are open. 2. Check the settings for flow and pressure. 3. Check if the tubing is pinched or crimped. 4. Check whether the pump housing locking device is engaged. 5. Check that the hose is guided correctly over the rollers.
Pump operates irregularly	<ol style="list-style-type: none"> 1. Reduce flow until pump operation stabilizes.
Inadequate pressure	<ol style="list-style-type: none"> 1. Choose higher pressure setting 2. Reduce the suction effect 3. Use high-flow cannulas
Distension liquid bloody or cloudy	<ol style="list-style-type: none"> 1. Increase the drain or suction effect of the Adapter Power System 2. Increase pressure 3. Increase the flow rate of the pump

5.2 Replacing the fuses

When changing the fuses, always ensure that only fuses with the same values are used. Replace the main fuse as follows:

1. Disconnect the device from the main power supply.
2. Apply a flat screwdriver to the slot at the bottom of the plug housing and gently pry it open.
3. Take the fuse holder out of the housing as far as possible and place it so that the fuses point upwards.
4. Gently lift the black lever between the fuses and pull the fuse bracket out.
5. Remove and replace the fuses.
6. Put the fuse bracket back in the holder with the contacts pointing outwards.
7. Put the fuse holder back in the plug housing.

5.3 Warranty

Arthrex GmbH assumes a warranty to the first purchaser for a twelve month period with regard to defects or failure of its instruments and electric devices. All defective products covered by the warranty are repaired or replaced free of charge by Arthrex at Arthrex's choice. The warranty does not cover any damage caused by unlawful use or improper handling of the particular product.

The warranty becomes invalid when Arthrex products are changed in any way or repairs are performed by anyone other than Arthrex.

Arthrex will gladly answer any questions referring to the quality, reliability and/or shelf life of this product.

6 TECHNICAL DATA

Control unit (AR-6450)

- Width: 36.8 cm (14.5 inches)
- Height: 12.7 cm (5 inches)
- Depth: 26.0 cm (10.25 inches)
- Weight: 6.6 kg (14.6 lbs.)
- Flow rate: 0 - 1600 ml/min
- Pressure: 0 - 255 mmHg
- Pressure control: continuous
- Operating mode: permanent
- Water protection: IPX1
- Main cable: 10 A/250 V
- Connector: CEE 7/7
- Jack: IEC 320/C13
- Power supply: 100-240 V, 50/60 Hz,
- Fuse: T2.0A
- Cleaning: Surface cleaning with mild detergent
- Sterilization: Surface disinfection with mild disinfectant

Remote control (AR-6306)

- Width: 63.5 mm (2.5 inches)
- Height: 95.3 mm (3.8 inches)
- Depth: 22.2 mm (0.9 inches)
- Weight: 0.23 kg (0.5 lbs.)
- Cable length: 3 m (9.8 feet)
- Cleaning: Surface cleaning with mild detergent

Ambient conditions for operation

- Temperature: 10° to 40°C (50° to 104°F)
- rel. humidity: 30% to 75%
- Air pressure: 700 hPa to 1060 hPa

Ambient conditions for storage (in shipping packaging)

- Temperature: -40° to +70°C (-40° to 158°F)
- rel. humidity: 0% to 100%, condensation

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© Copyright Arthrex GmbH, 2002. All rights reserved. LM6450 Rev. 1
Arthrex, Inc.
Arthrex S.A.S.
Lasswell Medical Co., Ltd., 405 Industrial Drive, Unit 21, Milton, Ontario • Canada L9T 5B1
Arthrex Ltd.
Dimco S.p.A., Via Degli Olimetti N. 56, 00060 Formello (RM) • Italy